


Placas Palmares


Perfil de la placa (mm)
Se muestran las placas izquierdas
I:izquierda
D:derecha, solo referencia



Placas Palmares

Placas Dorsales

Placas para minifragmentos



Alambres K acero inoxidable

Referencia	Sistema	∅	Descripción	Largo
A-5040.21		1.2 mm	trocar	150 mm
A-5040.41		1.6 mm	trocar	150 mm
A-5042.21		1.2 mm	lancet	150 mm
A-5042.41		1.6 mm	lancet	150 mm
A-5042.51		1.8 mm	lancet	150 mm
A-5040.00	2.2	0.8 mm	trocar	100 mm
A-5040.10	3.0	1.1 mm	trocar	100 mm
A-5042.00	2.2	0.8 mm	lancet	100 mm
A-5042.10	3.0	1.1 mm	lancet	100 mm
A-5043.00	2.2	0.8 mm	2 x trocar	100 mm
A-5043.10	3.0	1.1 mm	2 x trocar	200 mm

Tornillo Trilock Ø 2.5 HD7



Referencia	Largo
A-5750.08/1	8 mm
A-5750.10/1	10 mm
A-5750.12/1	12 mm
A-5750.14/1	14 mm
A-5750.16/1	16 mm
A-5750.18/1	18 mm
A-5750.20/1	20 mm
A-5750.22/1	22 mm
A-5750.24/1	24 mm
A-5750.26/1	26 mm
A-5750.28/1	28 mm
A-5750.30/1	30 mm

Tornillos en paquetes de 5

Tornillo Cortical Ø 2.5 HD7



Referencias	Largo
A-5700.08/1	8 mm
A-5700.10/1	10 mm
A-5700.12/1	12 mm
A-5700.14/1	14 mm
A-5700.16/1	16 mm
A-5700.18/1	18 mm
A-5700.20/1	20 mm
A-5700.22/1	22 mm
A-5700.24/1	24 mm
A-5700.26/1	26 mm
A-5700.28/1	28 mm
A-5700.30/1	30 mm

Tornillos en paquetes de 5

Brocas Stryker para tornillos de Ø 2.5 mm



A-3720
Ø 2.0 mm, tope 40 mm, largo 97 mm



A-3721
Ø 2.6 mm, tope 10 mm, largo 67 mm